



PROUDLY OFFERS

## SLM RPD

SELECTIVE LASER MELTED  
REMOVABLE PARTIAL DENTURE

- SLM- Selective Laser Melting  
Frameworks are 3D laser sintered directly to metal with unparalleled precision – no warpage!
- Fully CAD/CAM designed and manufactured frameworks
- 2 year full warranty
- 100% made in the USA with FDA approved materials (Chrome Cobalt Alloy)
- Precision fit and reliability
- High fatigue resistance

**PROFESSIONAL DENTAL COMPANY**

611 Omaha Street  
P.O. Box 1257  
Sioux City, Iowa 51102  
712-252-4034 • 800-831-0936  
Fax 712-252-3037

**JAEGER DENTAL COMPANY**

2218 West 2nd Street  
P.O. Box 4908  
Grand Island, NE 68802  
308-382-0235 800-289-3525  
Fax 308-382-2208

**BILLINGS DENTAL LABORATORY**

1407 N. 203rd Street  
Elkhorn, NE 68022-1771  
402-551-5954 • 800-747-5441  
Fax 402-551-0498